



Altogen Labs  
4020 South Industrial Dr  
Suite 130  
Austin TX 78744 USA  
Tel: (800) 658-7009

## **Cell Banking Services**

Our GLP-compliant cell line banking systems provide mammalian cell production, assurance that a uniform population of cells is preserved, cell integrity is maintained, and cells are readily available per customer request.

### **Facilities**

Cell line banks are maintained in a secure, controlled and monitored storage environment. Cell banking work is performed within certified Class II laminar flow biosafety cabinets with HEPA filtered air. Cells are grown only in pre-calibrated incubators with UV air and water protection.

### **Quality Control**

GLP-compliant process with established SOPs and associated cell line batch records. Cell lines are certified to be free of Mycoplasma or Bacterial contamination

### **Cell growth and Cryopreservation**

Altogen Labs provide cell expansion to manufacture a cell bank using standard T-flasks, roller bottles, shake flasks, and cell factories. Cryopreservation is performed in a monitored liquid nitrogen system. We can ship the cell bank (on dry ice) directly to the client or to a specified third party.

Altogen Labs can provide an immediate price quote (contact e-mail:[techserv@altogen.com](mailto:techserv@altogen.com)).

Please note that experimental details allow accurate quotes and timeline estimates.